

Corrigendum

Sl. No.	Nomenclature	Corrigendum Requested	Corrigendum
26	Ovum Pick up Needle	Single lumen separately; Double lumen separately; quote per piece	To be read as 26 A. Single lumen Ovum Pick up Needle Catheter Ovum Picks up Set (with valve) Tube 60- 100 cms Needle 30-35 cms. Gauge 17. Thin Echogenic tip Storage: Room temperature. Point 26 B Double Lumen Ovum Pick up Needle & Catheter Ovum Picks up Set (with valve) Tube 60- 100 cms Needle 30-35 cms. Gauge 17. Thin Echogenic tip Storage: Room temperature. Price to be quoted per piece for 26 A and B separately
27	Immature Ovum Pick up Needle	Single lumen separately; Double lumen separately; quote per piece	To be read as Ultrasound-guided transvaginal aspiration of immature oocytes for In Vitro Maturation (IVM) and fertility preservation procedures. 20-gauge needle. 30–35 cm needle length. Medical-grade stainless steel
28	Embryo Transfer Catheter	All 3 types separately; quote per piece	Point to be read as 28(A). Embryo Transfer Catheter . Soft (visible on USG). Outer sheath 17.3 (6.6F) cm. Inner 24 cm (2.8F) Storage : Room temperature. MEA tested
			28(B) Embryo Transfer Catheter Rigid (visible on USG). Length 23 cm. Inner catheter diameters

			1.52 mm (outer) and 0.76 mm (inner) Storage : Room temperature. MEA tested. Prices to be quoted per piece for each 28A and 28B separately.
35	Flexible Pipette	All tip sizes separately; quote per piece	Denuding cumulus mass from the oocytes prior to ICSI to check fertilization or for micromanipulation of oocytes embryos, blastocyst and blastomeres. Flexible polycarbonate pipette is designed for use with reusable pipetting handles. Supplied sterile, packing as – ten per vials.) Inner Diameter - 130/135 μm Inner Diameter - 140/150 μm Inner Diameter - 170/175 μm Inner Diameter - 225/250 μm Inner Diameter - 275/300 μm Prices to be quoted per vial for each size separately.
36	Handle for Flexible Pipette	Describe handle; quote per piece	Handle to be used with flexible pipettes for handling gametes and embryos. Aspiration volume = 0.25 microliter to 3 microliters. Accessory set consists of replacement part One small O ring, One large O ring, a spacer, and a collect.
38	Embryo Biopsy Pipettes/Needle	Day-5 biopsy pipette; quote per piece	Blastomere biopsy for micropipette 30°angled blastomere flat and polished (28-32 μm) for Day-5. Holding pipette. Tip 120 μm(OD)/ Tip 25μm(ID)/,35° bend. Storage: Room temperature.
52	Stripper Micropipette Handle	Quote per piece	Point to be deleted. To be read as point 36