

Corrigendum for 4K Ultra High Definition Arthroscopy System		
A)	Built in Wi-fi Router for wireless connectivity	Delete
B)	Powerful 300 Watt Xenon Lamp with emergency lamp facility	Delete
C)	Digital storage of still images, video sequences and audio files	Digital storage of still images, video sequences
	50-58" 4K monitor with Resolution of 3840 x 2160	32" 4k Monitor with resolution 3840X2160
	Camera Head should have 7 Years warranty against autoclaving	Camera Head should have 5 Years warranty against autoclaving
	Camera Head should be of Titanium Housing and Hermetically Sealed for Autoclaving	
D	Titanium housing kit	Titanium housing kit or equivalent
4)	Fiber Optic Light Cable	
	Must be Autoclavable	Must be Autoclavable/Soakable/ Plasma/ ETO
	3 Meter in length.	4-5 meters
5)	4K Arthroscopes 30°	
	4K Arthroscopes 30° should have diameter of 4.5 to 5 mm diameter with 152mm length	4K Arthroscopes 30° should have diameter of 4 to 5 mm diameter with 140-160 mm length
	4K Arthroscopes 30° should have diameter of 2.5 to 3 mm diameter with 152mm	4K Arthroscopes 30° should have diameter of 2 to 3 mm diameter with 120-140 mm length
	4K Arthroscopes 70° should have diameter of 4.5 to 5 mm diameter with 152mm length	4K Arthroscopes 70° should have diameter of 4 to 5 mm diameter with 140-160 mm length
7)	Arthroscopic Shaver System	
	Arthroscopic Shaver System consist of shaver console, hand piece and foot switch with option of connection Power Rasp....	Arthroscopic Shaver System consist of shaver console, hand piece and foot switch,
	Drill hand-piece	Drill console can be integrated with shaver or provided separately
	Max Speed 8000 RPM for forward and reverse movement and 3000 RPM for oscillating mode	Max Speed 12000 RPM for forward and reverse movement and 3000 RPM for oscillating mode
	Should set speed up/down	Delete
	OSC should have 3 modes(STD/AGG/Efficient)	Delete
		Handpiece Must be compatible with Small Joint and Large joint Blades or Separate hand piece should be provided.
9)	FLUSH button on the front panel	Optional.
	All quoted items should be of same manufacturer	All quoted items should preferable be of same manufacturer. Principal Company must take the warranty and CMC of the said third party products
	Consumables	
	Power rasp for SAD and Notchplasty	Delete
	Deflectable microfracture chondral pick...	Delete